


Instruction For Use B - Needle / Introducer



Step B-1

An epidural needle (8G) is inserted into the epidural space of the spine through Hiatus above the coccyx

Step B-2

Remove the Needle which is inside of epidural needle

Step B-3

Insert a flexible endoscope inside of V-catheter until it can get to the end of the V-catheter' Hub. Fill up V- catheter using a prefilled saline syringe.

Step B-4

insert the catheter into the targeted treatment point.

Step B-5

Check where the catheter is placed with C-Arm.

Inject the mixed medicine into the catheter through the injection hole of the body.


Ordering Information

Part No.	Description		Length	OD	EA/Kit
SWV30T	V-Catheter (Steerable Neurolysis Catheter)		530.0	Φ3.0	1EA
Part No.	Description		Length	Gauge	EA/Kit
SWVP-30V	Catheter Introducer Kit	Sheath	99.0	Φ4.0	1KIT
		Stylet	155.0	Φ3.0	
		Dilator	155.0	Φ4.0	
		Guide Wire	700.0	Φ0.85	
SWVP-8V	Needle/Introducer		152.0	Φ4.0	1EA

VEEELER

Steerable Neurolysis Catheter

Diagnosis and surgical procedure in the epidural space
Minimally invasive
Rapid recovery time



Steerable Neurolysis Catheter for back pain relief

The Steerable Neurolysis Catheter has been designed to allow the delivery of pain medication precisely within the lower spinal column. Appropriate for patients with chronic back pain not responsive to conservative treatment.



Indications

- A. Chronic Back Pain
- B. Post- Laminectomy Syndrome
- C. Spinal Stenosis
- D. Herniated Nucleus Pulposus

Feature

- Lightweight, ergonomic design
- Catheter outer diameter is 3mm
- Soft, atraumatic and radiopaque tip
- Radiopaque shaft
- Dual working channels
- Dual infusion ports with rubber buffers
- Direct visualization via endoscope
- 2-way steering
- Locking/unlocking button

Instruction For Use A – Catheter Introducer Kit

